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SAREMCO microcid etchant gel: Etching gel for selective enamel etching. The gel allows for visual in-process control of etching. Contains phosphoric acid. Corrosive!

James-2: light-curing, solvent-free single-component adhesive (primer and bond) for enamel and dentine.

A dental adhesive system, such as James-2 is used in combination with mhc to achieve a permanent connection between the restoration and tooth.

2. Composition:

mhc: inorganic fillers (barium glass and silica 76 % wt, 52 % by volume, particle size between 4–3000 nm), BisGMA, BisEMA, TEGDMA, catalysts, inhibitors, pigments

SAREMCO microcid etchant gel: water, phosphoric acid, gel former, colorant

James-2: hydroxyethylmethacrylate, urethanemethacrylate, polyalkenoate methacrylized, hydroxypropylmethacrylate, glycerinedimethacrylate, catalysts, inhibitors

3. Indication

1. Restoration of class I, II, III, IV and V cavities on anterior and posterior teeth
2. Extraorally tempered inlays and onlays (direct and indirect), veneers
3. Aesthetic corrections of interdental spaces, enamel hypoplasia, discolorations etc.
4. Blockings, fractures

4. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylate.

5. Side effects

In individual cases, contact allergies have been described in the case of products of a similar composition. To avoid pulp reactions, it is advisable to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

6. Interactions

Avoid underfilling materials which may hinder polymerisation owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol, fall under this category.

7. Processing stages

7.1. Tooth cleaning

Brush the tooth to be treated and the teeth neighbouring it with a fluoride-free toothpaste. Clean interdental spaces using strips and dental floss if necessary.

7.2. Shade selection

Under natural light, determine the shade using the SAREMCO shade guides (alternatively: VITA Lumin Vacuum shade guide).

7.3. Drying

Dry sufficiently, apply rubber dam.

7.4. Cavity preparation

Prepare the cavity in the usual way. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

7.5. Underfilling

In the case of preparations near pulp, it is advisable to protect the pulp using an underfilling (e.g. calcium hydroxide). Do not use zinc oxide eugenol!

7.6. Adhesion

7.6.1. Etch enamel with SAREMCO microcid etchant gel

Evenly apply SAREMCO microcid etchant gel to the dry enamel lamina. Allow to take effect for 30–60 seconds until you can see air bubbles forming. Rinse and dry thoroughly. Keep dry. Enamel that has been sufficiently etched has a matt and whitish appearance. Note: In dentine-bounded cavities, use the total etch technique if necessary. The dentine should be etched for no longer than 15 seconds and only dried briefly following rinsing.

7.6.2. Prime and bond with James-2

Apply James-2 to the cavity and rub in for 20 seconds until an evenly glossy film appears over the entire cavity. Light-cure for 30 seconds.

7.7. Restoration with mhc

Gradually adapt mhc up to a maximum layer thickness of 2 mm. Light-cure every layer for 40 seconds. When using matrices, after removing them expose filling to light again from lingual and buccal.

Warning: After each use, turn the syringe screw back by one turn and re-attach the screw cap. Protect light-curing products from strong sources of light. The specified exposure times refer to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm² and a wavelength of 400–500 nm. They apply to a maximum layer thickness of 2 mm and to all shades. The necessary exposure time may vary depending on the light source and its instructions for use. If in doubt, check the lamp's light output and the necessary exposure time before operation in vitro. During polymerisation, an inhibition layer, which must not be touched or removed if other composite layers are to be applied, forms on the surface.

Exposure time for increments of 2mm

Light power	≥ 500 mW/cm ²	≥ 800 mW/cm ²
mhc initial colours	20 sec.	20 sec.
mhc opaque colours	40 sec.	30 sec.
mhc all other colours	40 sec.	20 sec.

7.8. Finishing, polishing

Prepare the filling with 40 µ and 12 µ diamonds. Polish to a high-gloss using polishing brushes, polishing discs, strips or silicone polishers. The filling can be prepared and polished immediately after polymerisation.

8. Storage

Protect light-curing products from strong sources of light and heat! In the case of room temperatures above 28°C, it is advisable to store the products in the refrigerator. Prolonged temperatures above 28°C can shorten the shelf life of the product. Under these conditions, the shelf life is 48 months. The expiry date is specified on the packaging.

9. Batch number and expiry date

The batch number should be specified to identify products in the case of enquiries. Products should no longer be used once the expiry date has elapsed.

10. Precautionary measures

Close containers after each use using the right lid. Keep out of reach of children. When using tips, ensure before use that the cap has been removed.

11. Emergency measures

In the case of direct contact with the oral mucosa, rinse with water. In the case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

12. Hygiene

Use application instruments for one patient only. Clean and disinfect reusable instruments after use. Dispose products away from patients to avoid contamination.

13. Warranty

Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. We accept no liability for other damage, namely that caused by non-compliance with the instructions for use or other improper handling or unintended use of a product. Before using the products, the user should check whether they are suitable for the intended purpose. He alone assumes all of the risks associated with using the product and bears sole responsibility for any damage caused as a result of this.

14. Production and distribution

SAREMCO Dental AG
Gewerbestrasse 4
CH-9445 Rebstein/Switzerland
Tel: +41 71 775 80 90
Fax: +41 71 775 80 99

edited 2015-11 | D600158

Class IIa medical devices
Class III medical devices (Canada)

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