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1. Produktbeschreibung

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Instruction for use apt unibond

1. Product description

Light-curing, 1-component self-etching adhesive to create a permanent marginal-gap-free adhesion between the tooth structure (dentine, enamel) and the light-curing filling/filling material.

2. Composition

ethanol, methacrylated phosphoric salt, aromatic urethanmethacrylate, BisEMA, water, fillers, initiators

3. Indication

- Direct restorations with light-curing composite-based materials.
- Indirect restorations: at the use of light-curing composite cements to fix inlays, onlays, crowns and bridges.

4. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates.

5. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide). Note: apt unibond does not contain TEGDMA or HEMA.

6. Interactions

Avoid underfilling materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

7. Processing stages

7.1. Drying

Due to the following application of the filling material, a drainage of the working field is required. Avoid any contamination of the cavity with blood or saliva. A rubber dam is recommended.

7.2. Cavity preparation

Prepare the cavity as usual. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity. Avoid overdriving.

7.2.1. Etching of tooth substance (optional)

By applying apt unibond clinically sufficient bond strength values are achieved. By selective enamel etching or in total-etch technique even higher adhesion values can be achieved.

7.2.1.1. Total-etch technique (with SAREMCO cmf etc.)

Apply the etching gel and condition in correspondence to instruction for use of the product. SAREMCO recommends cmf etc.

7.2.1.2. Selective enamel etching

Apply the etching gel and condition in correspondence to instruction for use of the product. SAREMCO recommends cmf etc. Sufficiently etched enamel appears matt and white.

Note: Possibly use the total-etch-technique in cavities limited by dentine. In this case, the dentine should not be etched for more than 15 seconds and only dried briefly after rinsing.

7.2.2. Non-Etch

Do not overdry the dentine. A slightly moist cavity surface is optimal. Without preparation (e.g. cervical defect), thoroughly clean the surface and roughen the enamel margin with a finishing diamond.

7.3. Application

Apply apt unibond to enamel/dentine in a not too thin layer with a brush and massage into the tooth substance for 20 seconds.

Note: apt unibond is light-curing, so avoid too intense ambient light. Dim the surgical light during the application.

7.4. Drying

Dry the adhesive layer with an air blower for at least 5 seconds.

7.5. Polymerizing

Polymerize the bonding layer with blue light (halogen or LED lamp). The specified exposure time refers to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm² and a wavelength of 400 - 500 nm. The necessary exposure time may vary depending on the light source and its instructions for use. In case of doubt, check the lamp's light output and the necessary exposure time before operation in vitro. Bring the light curing device as close as possible to the surface to be polymerized. During polymerisation, an inhibition layer forms on the surface; this must not be touched or removed, since it is important for the bond with the subsequent composite layers.

Exposure time

light power	≥ 500 mW/cm ²	≥ 1'200 mW/cm ²
exposure time	20 sec.	10 sec.

Important: Extending the exposure times does not have a positive effect on the adhesion! When using dual-curing luting composites, light curing of the luting composite is absolutely necessary.

Note: Do not use apt unibond for opaque restorations or deep cavities where complete light-curing is not guaranteed.

7.6. Notes for fixing Maryland bridges

When fixing Maryland bridges with apt unibond, the enamel must be roughened orthched with phosphoric acid gel (e.g. SAREMCO cmf etc.) prior to use of the bond. Self-etching bonds are less effective on non-prepared enamel. Excessive amount of bond material onto unpolished and unbevelled enamel can cause edge staining. For optimal bonding results with apt unibond, the enamel should be prepared in advance in class III, IV, V and for diastema closures.

8. Storage

Close apt unibond bottle tightly after use. Store in an upright position to ensure the backflow of the liquid. Do not expose light-curing products to direct sunlight or operating light. apt unibond was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Constant temperatures above 28°C / 82°F can reduce the shelf-life of the product.

9. Batch number and expiry date

The batch number should be specified to identify products in the case of enquiries. Products should no longer be used once the expiry date has elapsed.

10. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes into contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice.

11. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

12. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination.

13. Warranty

Our liability is restricted to the quality of our products. In case of a product being of defective quality, only its value is replaced. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. It is the responsibility of a user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental.

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