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Instruction for use apt composite

1. Product description

Light-curing, highly filled, radio-opaque micro hybrid composite for anterior and posterior restorations with very high resistance to enzymatic degradation. DIN EN ISO 4049 compliant.

2. Composition

inorganic filler (barium glass and silica 74 % wt, 49 % by volume, particle size between 0.04 - 3.00 µm, median 0.7 µm), aromatic urethane methacrylate, BisEMA, silica, catalysts, inhibitors, pigments

3. Indication

1. Restoration of class I, II, III, IV and V cavities on anterior and posterior teeth
2. Extra orally tempered inlays and onlays (direct and indirect), veneers
3. Aesthetic corrections of interdental spaces, enamel hypoplasia, discolorations etc.
4. Blockings, fractures

4. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates.

5. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide). Note: apt composite does not contain TEGDMA or HEMA.

6. Interactions

Avoid underfilling materials which may hinder polymerisation owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

7. Processing stages

7.1. Tooth cleaning
Brush the tooth being treated and the teeth next to it with fluoride-free toothpaste. If necessary, clean interdental spaces using strips and dental floss.

7.2. Shade selection

SAREMCO shades relate to VITA standards.

7.3. Drying

Dry sufficiently, apply rubber dam.

7.4. Cavity preparation

Prepare the cavity as usual. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

7.5. Underfilling

In case of near-th pulp preparations, it is advisable to protect the pulp using an underfilling (e.g. calcium hydroxide). Do not use zinc oxide eugenol!

7.6. Adhesion

The adhesive has to be conditioned and applied in correspondence to instruction for use of the product. SAREMCO recommends cmf adhesive system, els duobond and els unibond.

7.7. Restoration with apt composite

Application of the tip: Remove tip cap. Push the tip into the corresponding recess of the injector and fix and click it into place. (SAREMCO recommends SAREMCO applicator). Apply material with constant pressure of the lever.

Gradually adapt apt composite up to a maximum layer thickness of 2 mm. For optimal adaptation titanium nitride coated instruments are recommended. Light-cure every layer for 40 seconds. When using matrices, after removing them expose filling to light again from lingual and buccal. Note: apt composite produces an extremely low shrinkage stress.

Warning: Protect light-curing products from strong sources of light.

The specified exposure times refer to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm² and a wavelength of 400 - 500 nm. They apply to a maximal layer thickness of 2 mm. The necessary exposure time may vary depending on the light source and its instructions for use. If in doubt, check the lamp's light output and the necessary exposure time before operation in vitro. During polymerisation, an inhibition layer forms on the surface; this must not be touched or removed, if other composite layers are going to be applied.

Exposure time for increments of 2 mm

light power	≥ 500 mW/cm ²	≥ 800 mW/cm ²
apt composite A1, A2, A3, A3.5/B4	40 sec.	20 sec.

7.8. Finishing, polishing

Prepare the filling with 40 µm and 12 µm diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers. The filling can be prepared and polished immediately after polymerisation.

8. Storage

Protect light-curing products from strong sources of light and heat! apt composite was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F freezing in the refrigerator is recommended. Do not freeze!

9. Batch number and expiry date

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

10. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. When using tips, make sure, that before usage the cap has been removed. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice. When polishing or removing composites, it is recommended to always use a water cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

11. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

12. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination. Sterilization of syringes or tips before first use is not required.

13. Warranty

Our liability is restricted to the quality of our products. In case of a product being of defective quality, only its value is replaced. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental.

14. Other notes for Europe

If the user and/or patient become aware of serious incidents connected with application the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides.

15. Production / distribution

SAREMCO Dental AG
Gewerbestrasse 4
CH-9445 Rebstein / Switzerland
Tel: +41 (0) 71 775 80 90
Fax: +41 (0) 71 775 80 99
info@saremco.ch
www.saremco.ch

Edited: 06/2021 | D600198

Class IIa medical devices
Class III medical device (Canada)

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