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### 1. Produktbeschreibung

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### 1. Product description

Flowable, light-curing, radio-opaque microhybrid lining-composite with extremely low shrinkage-stress and an increased curing depth (increment thickness 4 mm). DIN EN ISO 4049-compliant.

### 2. Composition

barium glass (53 % wt, 28 % by volume, particle size between 0.04 - 3.00 µm, median 0.7 µm), silanized, ytterbium trifluoride BisGMA, BisEMA, catalysts, inhibitors, additives

### 3. Indication

- Underfilling (filling base) and lining for class I and II cavities
- Direct filling therapy of small class I restorations in remaining molars and for class I & II restorations in primary molars (without the need for an additional occlusal composite material)
- Extended fissure sealing
- Core build-ups

### 4. Contra-indication

Opened pulp, pulpitis, and known allergy to methacrylates, camphorquinone and BHT.

### 5. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is advisable to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide). Note: els bulkfill does not contain TEGDMA or HEMA.

### 6. Interactions

Avoid underfilling materials which may hinder polymerisation owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

### 7. Processing stages

#### 7.1. Tooth cleaning

Brush the tooth being treated and the teeth next to it with fluoride-free toothpaste. If necessary, clean interdental spaces using strips and dental floss.

#### 7.2. Drying

Dry sufficiently, apply rubber dam.

#### 7.3. Cavity preparation

Prepare the cavity as usual. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

#### 7.4. Placing the matrices

Take the necessary precautions in areas where els bulkfill can flow off. The use of matrices that are wedged in the proximal area is an advantage. Minimal separation facilitates the approximal contact and the insertion of the matrice.

#### 7.5. Underfilling

In case of near-to-pulp preparations, it is advisable to protect the pulp using an underfilling (e.g. calcium hydroxide-cement). Do not use zinc oxide eugenol!

#### 7.6. Adhesion

The adhesive has to be conditioned and applied in correspondence to instruction for use of the product. SAREMCO recommends cmf adhesive system, els duobond, els unibond and apt unibond.

#### 7.7. Application

**Using the syringe:** Fix the application tube (included) onto the syringe and apply els bulkfill directly.

Starting with the lowest point, insert the material directly into the prepared cavity and fill the cavity from bottom to top. Leave the tip of the metal tube inserted in the cavity. To prevent bubbles from forming, make sure that the material flows slowly and evenly. If necessary, remove bubbles and excess material before polymerisation.

When using as a basis for filling, the distance to the occlusal edge of the cavity must be at least 2 mm. els bulkfill can be polymerised in increments of up to 4 mm thick. Then cover with a methacrylate-based universal or posterior tooth composite. SAREMCO recommends els composite. Please note the specific instructions for use for processing and finishing.

**Note:** After each use, remove the application tube and put the cap back in place. Light-curing products should be kept away from strong light sources.

### 7.8. Light-curing

Use standard polymerisation devices for light-curing. Expose every area of the filling sufficiently. Additional exposure from the approximal space (where accessible) and the oral and vestibular spaces is also recommended. If metal matrices have been used, they must be removed beforehand. Bring the light window of the lamp as close as possible to the surface of the filling (otherwise the curing quality will be poorer). Insufficient curing can cause discoloration of the material and loss of strength of the composite filling and can lead to complaints.

The exposure time is 20 seconds (for halogen lamps of 700 mW/cm<sup>2</sup>). The exposure time may vary depending on the light source and its instructions for use.

### 8. Storage

Protect light-curing products from strong sources of light and heat! els bulkfill was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Bring to room temperature before use.

### 9. Batch number and expiry date

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

### 10. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice. When polishing or removing composites, it is recommended to always use a water cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

### 11. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

### 12. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination.

### 13. Warranty

Our liability is restricted to the quality of our products. In case of a product being of defective quality, only its value is replaced. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental.

### 14. Other notes for Europe

If the user and/or patient become aware of serious incidents connected with application of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides.

### 15. Production / distribution

SAREMCO Dental AG  
Gewerbstrasse 4  
CH-9445 Rebstein / Switzerland  
Tel: +41 (0) 71 775 80 90  
Fax: +41 (0) 71 775 80 99  
info@saremc.ch  
www.saremc.ch

Edited: 01/2021 | D600154

Class III medical devices  
Class III medical device (Canada)

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