

Instruction for Use I saremco print - CROWNTEC I US-version

**1. Product description**  
**saremco print CROWNTEC** is a light-curing, flowable polymer based on methacrylic acid ester for production of 3D-printed permanent crowns, inlays, onlays, veneers, temporary crowns & bridges, and artificial teeth (i.e., complete or partial dentures).

**saremco print** products are part of an overall concept of 3D printable resin-based materials and may only be used in combination with the specified printers and recommended equipment and in compliance with the manufacturer's instructions.

**Note** – The use of noncompliant devices might impair the function of the restoration. Sole responsibility for correct application is assumed by the user and is beyond control of SAREMCO Dental AG. SAREMCO Dental AG does not assume any responsibility and liability for damages caused by misuse.

**2. Composition**  
Enrichment products of 4,4'-isopropylidenediphenol, ethoxylated and 2-methylprop-2-enic acid, silanized dental glass, Pyrogenic silica, initiators. Total content of inorganic fillers (particle size 0.1 µm) is 30 - 50 % by mass.

**3. Intended Use**  
**saremco print** products provide light-curing 3D printable resin-based materials for the correction or reconstruction of functionally compromised natural dentition (e.g., missing teeth or deficient teeth) by manufacturing of customized 3D-printed dental prostheses.

**saremco print CROWNTEC** is to be used with 3D-printers from Nextdent, ASIGA, Rapid Shape, SprintRay, Prozen or Accuretta for the following applications (see section 8 herein):  
• Production of permanent crowns, inlays, onlays and veneers  
• Production of temporary crowns and bridges, inlays, onlays and veneers  
• Production of artificial teeth for subsequent insertion into a denture base

**4. Indications for Use**  
**saremco print CROWNTEC** is a light-curing 3D-printed material intended as an indirect restorative for both anterior and posterior restorations, including occlusal surfaces. The **CROWNTEC** material is used for fabricating permanent restorations such as inlays, onlays, veneers and full crown restorations. **saremco print CROWNTEC** can also be used for the fabrication of artificial teeth and temporary crowns & bridges.

**5. Contraindications**  
**saremco print CROWNTEC** is contraindicated for the following applications:  
• Maryland bridge, Inlay bridge  
• All forms of cantilever bridges  
• Bruxism of the patient  
• More than one pontic  
• Do not use the product in case of a known allergy to one or more ingredients.  
• In case of doubt, clarify and exclude a possible allergy with the help of a specific allergy test before using **saremco print CROWNTEC**.  
**saremco print CROWNTEC** must not be used for any other purposes than those specified in the "Indication" section. Any deviation from the instructions for use may have negative effects on the chemical and physical quality of the restorations produced from **saremco print CROWNTEC**.

**6. Interactions**

**Note** - known.

**7. Material Properties**

Color**	A1, A2, A3, B1, sw	Flexural strength**	a 120 MPa Average ± 135 MPa
Density**	ca. 1.4-1.5 g/cm <sup>3</sup>	Layer thickness when printing	50 µm
Viscosity*	2300-6000 mPa·s	Wavelength 3D-printer	365 or 405 nm

\*Applies to liquid resin

\*\*Applies to cured plastic printed with a 3D-Printer

**8. Requirements**

**Printers**  
• Nextdent S100 Figure 4 (405 nm)  
• Asiga MAX UV & PRO 4K (385 nm)  
• Rapid Shape D20 II, D30 II, D40 II, D10+, D20+ -cartridge, D20+, D30+ & D40+ (385 nm)  
• SprintRay K55 & K95 (405 nm)  
• Prozen Sonic XL 4K & Sonic 4K (405 nm)  
• Accuretta SOL & DENTIQ (405 nm)

**Software**

• Autodesk Netfabb  
• Composer  
• 3D Sprint  
• Rayware  
• DS Slicer  
• Alpha 3D

**Post Curing Unit**


• Ctoflash G171 (NIK-optik)  
• Signum HiLite Power (Kulzer)  
• Ctoflash Pro (Nextdent)  
• Prozen Cure V2  
• Cune (Accuretta)

**9. Processing Steps**

The following instructions have to be observed during tooth preparation:

Make sure to avoid tangential, spring edge or lip preparations as they are contraindicated with printed restorations. Therefore, exercise special care when using instruments with a round tip and do not introduce them any further than up to half their diameter at maximum. Please note that tangential preparations are technically unfeasible and would result in too thin, i.e., unstable and over-contoured, crown margins.

The following instructions apply to the model modeled on the computer:  
**Minimum Wall Thickness** – The following illustration shows the specified minimum wall thicknesses for the respective indication. The wall thickness must not be undercut even after manual grinding.



The following applies to temporary bridges: connector area at least 16 mm<sup>2</sup>. The connector area should be as large as possible. For physical stability, the height of the connector is more important than the width. The connector should be only doubling the strength, while doubling the height results in eight times the strength. Oval connector area are three-fold recommended.

Make sure that enough supports are generated. It is recommended to place the supports on the occlusal surface.

**9.1. Generating Printing File**  
Generate the printing file of the desired restoration by using appropriate software (Autodesk Netfabb, Composer, 3D Sprint, Rayware, DS Slicer or Alpha 3D) and deliver it suitable to the printer. Please observe the corresponding instruction for use of software and printer. Select the build style / IN file / material file for **saremco print CROWNTEC** in the printer software. Make sure that all software is up to date.  
**Important Note** – Commercially available artificial teeth may be subject of copyright law. When using an STL file of such third party, copyright laws must be considered.

**9.1.1 Nextdent Printer (S100 Figure 4) and Software**

**I. Hardware**  
Please refer to the printer's manufacturer's manual for this information.  
See the applicable user guides (<http://infocenter.3dsystems.com/nextdentS100-user-guide>)

**II. Nextdent printer software – 3D Sprint**  
Please refer to the printer's manufacturer's manual for this information.  
See the applicable user guides (<https://support.3dsystems.com/articles/3D-Sprint>)

**III. Printing parameters**  
Printing parameters are automatically loaded into 3D Sprint.  
a. Support parameters are automatically loaded into 3D Sprint.  
b. Slice thickness: 50 µm  
c. Optimal orientation: 0 degree tilted orientation  
d. Support parameters are automatically generated in 3D Sprint.  
**Note** – The occlusal side must face the build platform.

**IV. Environmental Conditions**  
a. Temperature of 3D-printing should be kept at 18 - 28°C (64.4 - 82.4°F)  
b. Humidity: 20 - 70 %

**9.1.2. ASIGA Printer (MAX UV & PRO 4K) and Software**

**I. Hardware**  
Please refer to the printer's manufacturer's manual for this information.

**II. Asiga printer software – Composer**  
Please refer to the printer's manufacturer's manual for this information.

**III. Printing parameters**  
Downloaded the required parameter set from the ASIGA database. A working temperature of 35°C / 95°F must be maintained.  
a. Slice thickness: 50 µm  
b. Optimal orientation: 0 degree tilted orientation  
c. Support point size: varies based on support type chosen.  
d. Support density: perimeter of the restoration and occlusal region.  
**Note** – The occlusal side must face the build platform.

**IV. Environmental Conditions**  
a. MAX UV print temperature: 35 ± 3 °C // 95 ± 3 °F  
b. Humidity: 20 - 80 %

**9.1.3 Rapid Shape Printer (D20 II, D30 II, D90 II) and software**

**I. Hardware**  
Please refer to the printer's manufacturer's manual for this information.

**II. Rapid Shape printer software – Autodesk netfabb**  
Please refer to the printer's manufacturer's manual for this information.

**III. Printing parameters**  
a. Slice thickness: 50 µm  
b. Optimal orientation: 0 degree tilted orientation  
c. Support point size: varies based on support type chosen.  
d. Support density: perimeter of the restoration and occlusal region.  
**Note** – The occlusal side must face the build platform.

**IV. Environmental Conditions**  
a. Print temperature: Room temperature; Please refer to the printer's manufacturer's manual for this information.  
b. Humidity: 20 - 80 %

**9.1.4 SprintRay Printer (K55 and K95) and software**

**I. Hardware**  
Please refer to the printer's manufacturer's manuals for this information.

**II. SprintRay printer software – Rayware**  
Please refer to the printer's manufacturer's manual for this information.

**III. Printing parameters**  
a. Slice thickness: 50 µm  
b. Optimal orientation: 0 degree tilted orientation  
c. Support point size: varies based on support type chosen.  
d. Support density: perimeter of the restoration and occlusal region.  
**Note** – The occlusal side must face the build platform.

**IV. Environmental Conditions**  
a. Print temperature: Room temperature; Please refer to the printer's manufacturer's manual for this information.  
b. Humidity: 20 - 80 %

**9.1.5. Prozen Printer (Sonic XL 4K & Sonic 4K) and software**

**I. Hardware**  
Please refer to the printer's manufacturer's manuals for this information.

**II. Prozen printer software – DS Slicer**  
Please refer to the printer's manufacturer's manual for this information.

**III. Printing parameters**  
Download the required parameter set from the Prozen database.  
a. Slice thickness: 50 µm  
b. Optimal orientation: 0 degree tilted orientation  
c. Support point size: varies based on support type chosen.  
d. Support density: perimeter of the restoration and occlusal region.  
**Note** – The occlusal side must face the build platform.

**IV. Environmental Conditions**  
a. Print temperature: Room temperature 25 ± 3 °C // 77 ± 3 °F; Please refer to the printer's manufacturer's manual for this information.  
b. Humidity: 20 - 80 %

**9.1.6. Accuretta Printer (SOL & DENTIQ) and Software**

**I. Hardware**  
Please refer to the printer's manufacturer's manuals for this information.

**II. Accuretta printer software – Alpha 3D**  
Please refer to the printer's manufacturer's manual for this information.  
Please refer to the Alpha3D video playlist information (<https://www.youtube.com/watch?v=frYfMwWmSc&list=PLYr-3ONrGc5Sm5dSdGvDnUvDq2d>)

**III. Printing parameters**  
Download the required parameter set from the Accuretta database.  
a. Slice thickness: 70 µm  
b. Optimal orientation: 0 degree tilted orientation  
c. Support point size: varies based on support type chosen.  
d. Support density: perimeter of the restoration and occlusal region.  
**Note** – The occlusal side must face the build platform.

**IV. Environmental Conditions**  
a. Print temperature: Room temperature 20 - 28 °C // 36 - 50 °F; Please refer to the printer's manufacturer's manual for this information.  
b. Humidity: 20 - 80 %

**9.2. Printing**  
Work as clean as possible, as dirty reservoirs or machines can cause deformation/discoloration and therefore failure of the printed objects.  
Briefly shake the liquid material and pour it into the reservoir of the 3D-printing machine. Start the printing process by following the instruction for use of the printer.  
**Caution** – Any unauthorized changes to the process equipment, parameters, or software may result in a device that is out of specifications. This is explicitly not recommended and is the responsibility of the user. In case of questions the user should contact the manufacturer for a list of validated software and process hardware.

**9.3. Cleaning**  
To achieve the desired material properties and biocompatibility, post-curing of the completely dried and cleaned printed objects is necessary. For final polymerization place the printed jobs in a UV-light box.  
**Note** – time of curing depends greatly on type of lamps / lightbox used. The final properties and the final color depend on the post-curing process. Post-curing is a UV-light treatment to ensure that **saremco print** materials obtain full polymer conversion, the residual monomer is reduced to a minimum and the highest mechanical properties are achieved.  
This procedure is a necessary step to attain a biocompatible end-product.  
It is suggested to use the polymerization unit "Signum HiLite Power" from (2 x 180 sec., turn around after 180 sec.) or the UVFlash device "Ctoflash G171" from NIK-optik (2 x 2000 flashes, turn around after 2000 flashes, UV-lbow: bowl and nitrogen). In general, lightboxes for light-curing veneering materials can be used that cover a wavelength range of 300 - 500 nm.  
Lightboxes with integrated flashlight allow shorter exposure time compared to conventional lamps. Always follow the respective instruction for use of the polymerization unit.  
Blast the surface of the printed jobs with blast polishing material carefully (e.g., Perlablast micro BEGO). Afterwards remove the support structures by using a cut-off-wheel or a cutter.

**Warning** – Protect light-curing products from strong light sources.

**9.4. Finishing the printed jobs**

To achieve the desired material properties and biocompatibility, post-curing of the completely dried and cleaned printed objects is necessary. For final polymerization place the printed jobs in a UV-light box.  
**Note** – time of curing depends greatly on type of lamps / lightbox used. The final properties and the final color depend on the post-curing process. Post-curing is a UV-light treatment to ensure that **saremco print** materials obtain full polymer conversion, the residual monomer is reduced to a minimum and the highest mechanical properties are achieved.  
This procedure is a necessary step to attain a biocompatible end-product.  
It is suggested to use the polymerization unit "Signum HiLite Power" from (2 x 180 sec., turn around after 180 sec.) or the UVFlash device "Ctoflash G171" from NIK-optik (2 x 2000 flashes, turn around after 2000 flashes, UV-lbow: bowl and nitrogen). In general, lightboxes for light-curing veneering materials can be used that cover a wavelength range of 300 - 500 nm.  
Lightboxes with integrated flashlight allow shorter exposure time compared to conventional lamps. Always follow the respective instruction for use of the polymerization unit.  
Blast the surface of the printed jobs with blast polishing material carefully (e.g., Perlablast micro BEGO). Afterwards remove the support structures by using a cut-off-wheel or a cutter.

**Caution** – Wearing nitrile gloves, safety glasses and dust mask is advised during this finishing process.

**Recommended light curing equipment (Post curing process)**

Manufacturer & Model	Ctoflash G171 (NIK-optik)	Signum HiLite Power (Kulzer)
Supply voltage	100, 117, 230 V	100, 115, 230 V
Lamp power	2 x 100 W	200 W
Light intensity	Not adjustable	Not adjustable
Wavelength	280 - 700nm	320 - 540 nm
Curing time	2 x 2000 flashes, turn around after 2000 flashes, "UV-lbow bowl" and nitrogen	2 x 180 sec., turn around after 180 sec

Manufacturer & Model	Nextdent LC-3D Printbox	Prozen Cure V2
Supply voltage	AC110-240V, 50/60 Hz	100-240V, 50-60Hz
Lamp power	Max. 264 W	60 W in total
Light intensity	Not adjustable	Not adjustable
Wavelength	Full light spectrum (300-550 nm)	365, 385, 405 nm
Curing time	30 minutes Optional: 2 min. boiling water bath after post curing	(1) 2 side x 5 min. per side curing (turn print in-between) (2) 2 min. boiling water bath after post curing

Manufacturer & Model	Cune (Accuretta)
Supply voltage	AC110-240V 50/60 Hz
Lamp power	70 W
Light intensity	11225 mW/cm <sup>2</sup>
Wavelength	365, 385, 405 nm
Curing time	T = 2 x 3 minutes P = 16 B = 10 B = ON 2 min. boiling water bath after post curing

**9.5. Fastening**

**9.5.1. Fastening the definitive crowns, inlays, onlays and veneers**  
In case of definitive single crowns, the inside of the crowns should be roughened with a sandblast (Al<sub>2</sub>O<sub>3</sub> 110 µm). Then, as usual, fix it definitively with a composite cement material. Zinc-phosphate cements as well as glass-ionomer-cements are only of limited suitability, due to their opacity. The fastening composites Panavia V5 [Kuraray] and Variolink [3M] are recommended.

**9.5.2. Fastening the temporary crowns and bridges, inlays, onlays and veneers**  
Fasten the finished temporary prosthesis with commercially available provisional cements.

**9.5.3. Connecting of artificial teeth and prosthesis**  
**a) Inserting the printed artificial teeth in a printed, prefabricated denture base.**  
Roughen the base surface of the printed artificial teeth for example by sandblasting (Al<sub>2</sub>O<sub>3</sub> 110 µm), apply a primer and a fixing material, insert it completely according to the natural shape and polymerise.

Alternatively, **saremco print CROWNTEC** can also be used directly as fixing material. Therefore place a small amount of material with a brush on the roughened teeth-surface of the artificial tooth, place it into the prosthesis, eliminate any excess material and light cure it from all sides for at least 20 seconds. The polymerisation light should have at least a light output of 600 mW/cm<sup>2</sup>.  
**b) Classical pouring method with cold cure resin.**  
After roughening the teeth, follow classical finishing procedures.

**9.6. Finishing, polishing**  
Prepare the restoration with 40 µ and 12 µ diamond burrs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishes.

**9.7. Additional advice**  
Remove the container from the printer and filter the resin through a fine 190 µm particle strainer, if:  
• printed material is not completely cured according to the - particulates of polymerized residues are visible in the container or stick to the bottom.  
Discard and replace the **CROWNTEC** material with a new batch, if discoloration, evident gelation, or polymerization is observed after filtering.

Do not mix different batches of **CROWNTEC** material.

**10. Storage**  
Protect this product from strong light and heat sources. The recommended storage temperature is between 4°C and 28°C / 39°F and 82°F. Close the package after each use.

**11. Batch number and expiry date**  
The batch number is used to identify the product in case of queries. Do not use this product after the expiration date.

**12. Precautionary measures**  
For medical use only. Keep out of reach of children.  
• The use of nitrile gloves while working with **saremco print CROWNTEC** is recommended until post-curing. Commercially available nitrile gloves do not provide effective protection against the sensitizing effect of methacrylates. If the product comes into contact with the glove, remove the glove and dispose of it, wash your hands immediately with soap and water and put on a new glove.  
• In case of an allergic reaction, consult a doctor.  
• When polishing or removing composites, it is recommended to always use a water-cooling system and a good extraction system, to ventilate the dental laboratory and to wear masks with high particle filtration efficiency for small particle sizes.

**Warnings**  
Hazardous components: ethoxylated Bisphenol A dimethacrylate  
H315 causes skin irritation / H317 may cause an allergic skin reaction / H319 causes serious eye irritation / H335 may cause respiratory irritation

**Precautionary Statements**  
P261 avoid breathing dust/fume/gas/mist/vapors/spray  
P274 wash with water and soap thoroughly after handling  
P271 use only outdoors or in a well-ventilated area  
P272 contaminated work clothing should not be allowed out of the workplace  
P273 avoid release into the environment  
P303+P361+P531 IF ON SKIN: wash with plenty of water  
P333+P313 if skin irritation or rash occurs: get medical advice/attention  
P362+P354 take off contaminated clothing and wash it before reuse  
P304+P340 IF INHALED: remove person to fresh air and keep comfortable for breathing  
P312 call a POISON CENTER/doctor if you feel unwell  
P305+P351+P338 IF IN EYES: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing  
P337+P313 if eye irritation persists: get medical advice/attention  
P403+P233 store in a well-ventilated place. Keep container tightly closed  
P501 dispose of contents/container to residual waste as per local and national regulations

**13. Emergency Measures**  
• In case of direct contact of the uncured material with the oral mucosa, rinse with water.  
• In case of contact with eyes, rinse thoroughly with water and consult an eye specialist.  
• In case of swallowing the restoration, consult a clinician.  
• In case of breakage or fracture of the restoration, consult a dentist.  
• In case of bleeding caused by the dental restoration, consult a clinician.  
• In case of an infection beneath dental restoration caused by the product, consult a dentist.

**14. Hygiene**  
Restorations made of **saremco print CROWNTEC** should not be cleaned with chemical products. Cleaning with water is sufficient. The finished restorations can - if necessary - be disinfected with an ethanol solution.

**15. Warranty**  
The product is developed for use in dentistry and must be processed in accordance with the instructions for use. For further details, please refer to the instructions for use. For the use of the product, it is necessary to comply with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

Scope of delivery	Contents	Packaging	REF
CROWNTEC_A1	300 g	Bottle	8063
CROWNTEC_A2	300 g	Bottle	8052
CROWNTEC_A3	500 g	Bottle	8051
CROWNTEC_B1	500 g	Bottle	8065
CROWNTEC_SW	500 g	Bottle	8066

**17. Production / distribution**  
SAREMCO Dental AG  
Gewerbestraße 4  
CH-9445 Rebasen / Switzerland  
Tel: +41 (0) 71 775 80 90  
Fax: +41 (0) 71 775 80 91  
info@saremco.ch  
www.saremco.ch

Edited: 02/2022 | D600248  
Saremco AG is a medical device  
US: Class II medical device

**Glossary**

	Manufacturer		Protect from Sunlight
	Batch Code		Temperature Limitation
	Item Number		Expiry Date
	CE Marking of Conformity		Use only by Professionals

Prescription Only

Rx Only