Gebrauchsanweisung els seal

1. Produktbeschreibung

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Light-curing, flowing fissure sealing resin. Compliant with DIN EN ISO 6874.

2. Composition

BisEMA, dental glas (silanized), silicium dioxid, catalysts, inhibitors

3. Intended Use

Saremco restoratives are intended for the direct reconstruction or correction of functionally compromised natural dentition (e.g., deficient teeth).

4. Indication

1. Preventive sealing of untreated or extended fissures or pits on milk teeth and permanent teeth
2. Sealing/veneering of damaged enamel surfaces
3. Anchoring of orthodontic regulatory systems and covering of predilection sites of caries in orthodontic treatments
4. Sealing of the marginal gaps in restorations made of composite, amalgam, ceramic

5. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates.

6. Side effects

In individual cases, contact allergies have been described to products of a similar composition.

Note: *els seal* does not contain TEGDMA, HEMA or BisGMA.

7. Interactions

None known.

8. Processing stages

8.1. Fissure cleaning

Brush the tooth to be treated and its neighboring teeth with a fluoride-free tooth cleaning paste. Thoroughly clean and rinse fissures with an airflow powder jet device with sodium bicarbonate or with a brush and pumice stone.

8.2. Fissure enlargement

If necessary, slightly enlarge fissures with a suitable instrument.

8.3. Enamel etching

Apply the etching gel and condition in correspondence to instruction for use of the product. SAREMCO recommends *cmf etch*. Rinse thoroughly, dry and keep dry. Sufficiently etched enamel appears matt and white.

8.4. Seal and cure

Remove syringe cap and screw on luer-lock cannula. Apply material with constant pressure on the syringe plunger. Do not retract the plunger after application (to prevent a possible slumping of the material). A retraction of the plunger causes air to be sucked into the material, thereby deteriorating the material properties. Light cure *els seal* for 20 seconds (halogen lamps with 700 mW/cm²) or 10 seconds (LED lamps with 1100 mW/cm²) per section. The necessary exposure time may vary depending on the light source and its instructions for use. In case of doubt, check the lamp's light output and the necessary exposure time before operation in vitro.

Warning: After each usage, the application tip must be removed, and the syringe tightly resealed with the original cap. Protect light curing products from strong sources of light.

8.5. Finishing, polishing

If necessary, re-seal defective areas and remove excess material. Possibly polish the seal with the aid of silicone polishers. It is recommended to fluorideate the treated teeth.

9. Storage

Protect light-curing products from strong sources of light and heat! *els seal* was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Bring to room temperature before use.

10. Batch number and expiry date

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

11. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. For dental use only. Commercially available medical gloves do not provide protection against the sensitization effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice. When polishing or removing composites, it is recommended to always use a water cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

12. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

13. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination. Sterilization of syringes or tips before first use is not required.

14. Warranty

The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

15. Other notes for Europe

If the user and/or patient become aware of serious incidents connected with application the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides.

16. Production / distribution

SAREMCO Dental AG
Gewerbstrasse 4
CH-9445 Rebstein / Switzerland
Tel: +41 (0) 71 775 80 90
Fax: +41 (0) 71 775 80 99
info@saremco.ch
www.saremco.ch

Edited 03/2022 | D600054

Class IIa medical devices
Class III medical devices (Canada)

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