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2. Composition

Inorganic fillers (barium glass and silica 76 % wt, 52 % by volume, particle size between 0.04 - 3.00 µm, median 0.7 µm), BisGMA, BisEMA, TEGDMA, catalysts, inhibitors, pigments

3. Intended Use

Saremco restoratives are intended for the direct reconstruction or correction of functionally compromised natural dentition (e.g., deficient teeth).

4. Indication

1. Restoration of class I, II, III, IV and V cavities on anterior and posterior teeth
2. Extra orally tempered inlays and onlays (direct and indirect), veneers
3. Aesthetic corrections of interdental spaces, enamel hypoplasia, discolorations
4. Blockings, fractures

5. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates.

6. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

7. Interactions

Avoid underfilling materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

8. Processing stages

8.1. Tooth cleaning

Brush the tooth being treated and the teeth next to it with fluoride-free toothpaste. If necessary, clean interdental spaces using strips and dental floss.

8.2. Shade selection

SAREMCO shades relate to VITA standards.

8.3. Drying

Dry sufficiently, apply rubber dam.

8.4. Cavity preparation

Prepare the cavity as usual. Undercuts and beveled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

8.5. Underfilling

In case of near-to-pulp preparations, it is advisable to protect the pulp using an underfilling (e.g. calcium hydroxide). Do not use zinc oxide eugenol!

8.6. Adhesion

The adhesive has to be conditioned and applied in correspondence to instruction for use of the product. SAREMCO recommends **James-2, els unbond** and **cmf adhesive system**.

8.7. Restoration with mhc composite

After each usage, turn the syringe screw back by one turn and re-attach the screw cap. Gradually adapt **mhc composite** up to a maximum layer thickness of 2 mm. When using matrices, remove them and expose filling to light again from lingual and buccal.

Warning: Protect light-curing products from strong sources of light.

The specified exposure time refers to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm² and a wavelength of 400 - 500 nm. They apply to a maximal layer thickness of 2 mm. The necessary exposure time may vary depending on the light source and its instructions for use. In case of doubt, check the lamp's light output and the necessary exposure time before operation in vitro. During polymerization an inhibition layer forms on the surface; this must not be touched or removed if other composite layers are going to be applied.

Exposure time for increments of 2 mm

light power	≥ 500 mW/cm ²	≥ 800 mW/cm ²
mhc opaque color	40 sec	30 sec
mhc all other colors	40 sec	20 sec

8.8. Finishing, polishing

Prepare the filling with 40 µm and 12 µm diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers. The filling can be prepared and polished immediately after polymerization.

9. Storage

Protect light-curing products from strong sources of light and heat! **mhc composite** was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Bring to room temperature before use.

10. Batch number and expiry date

The batch number should be specified to identify products in the case of enquiries. Products should no longer be used once the expiry date has elapsed.

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SAREMCO Dental AG