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1. Produktbeschreibung

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Instruction for use els duobond

1. Product description

Dual-curing two-component self-etching bond to create a permanent marginal-gap-free adhesion between the tooth structure (dentine, enamel) and the light-, self- or dual-curing filling/fixing material.

2. Composition

els duobond base: BisEMA, methacrylated phosphoric salt, initiators
els duobond catalyst: acetone, water, initiators

3. Intended Use

Saremco dental bonding agents promote adhesion of direct or indirect prosthesis for reconstruction or correction of functionally compromised natural dentition (e.g., deficient teeth).

4. Indication

- Direct self- or dual-curing composite restorations and core build-ups
- Direct restorations with light-curing composite-based materials
- Fixing root pins with dual- or self-curing composite cements
- Indirect restorations: at the use of dual-curing or self-curing composite cements to fix inlays, onlays, crowns and bridges.

5. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates.

6. Patient target group:

els duobond is suitable for use in all patients without any age or gender restrictions.

7. User

els duobond should only be used by a professionally trained dental practitioner.

8. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

Note: els duobond does not contain TEGDMA or HEMA.

9. Interactions

Avoid underfilling materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol, belong to this category.

10. Processing stages

10.1. Drying

Due to the following application of the filling material, a drainage of the working field is required. Avoid any contamination of the cavity with blood or saliva. A rubber dam is recommended.

10.2. Cavity preparation

Prepare the cavity as usual. Undercuts and beveled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity. Avoid overdrilling.

10.3. Mixing

els duobond is dual-curing and therefore must be processed immediately after mixing. Thoroughly mix 1 drop of the base and 1 drop of the catalyst on a mixing tray with a brush for around 5 seconds (forms a self-etching adhesive). After use, rub the adhesive seal with alcohol to avoid contamination.

10.4. Application

Apply els duobond to enamel/dentine in a not too thin layer with a brush and massage into the tooth substance for 20 seconds. **Note:** els duobond is light curing, so avoid too intense ambient light. Dim the surgical light during the application.

10.5. Drying

Dry the adhesive layer with an air blower for at least 5 seconds.

10.6. Polymerizing

10.6.1. Direct restorations

When using light-, self- or dual-curing filling material polymerize the bonding layer for 10 seconds with blue light (halogen or LED light).

10.6.2. Indirect restorations

Light-curing luting composites: When using light-curing luting composites polymerize the bonding layer for 10 seconds with light (halogen or LED light). **Self- or dual-curing luting composites:** When using dual- or self-curing luting materials, els duobond creates good adhesion even without light-curing. Optionally light-curing can be additionally performed for 10 seconds with blue light (halogen or LED lamp).

10.7. Fixing root pins

10.7.1. Bonding

To fix root pins, the mixed els duobond (base and catalyst) is massaged into the root canal and onto the occlusal surface of the stump with an applicator for 20 seconds. Disperse solvent with oil free air for 5 seconds. Remove any excess material by dabbing with an absorbent paper strip. Do not light-cure the bonding layer!

10.7.2. Inserting

Insert the root pin prepared according to the manufacturer's instructions into the root canal using a suitable fixing material in accordance with the instructions for use and to achieve slight excess cement. SAREMCO recommends els cem.

10.7.3. Light-curing

Light cure for at least 40 seconds to fix the root pin and to harden the occlusal applied els duobond.

10.8. Notes for fixing Maryland bridges

When fixing Maryland bridges with els duobond, the enamel must be roughened or etched with phosphoric acid gel (e.g. SAREMCO cmf etch) prior to use of the bond. Self-etching bonds are less effective on non-prepared enamel. Excessive amount of bond material onto unpolished and unbeveled enamel can cause edge staining. For optimal bonding results with els duobond, the enamel should be prepared in advance in class III, IV, V and for distal areas closures.

11. Storage

Close els duobond bottles tightly after use. Store in an upright position to ensure the backflow of the liquid. Do not expose light-curing products to direct sunlight or operating light. els duobond was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store in a refrigerator at 4°C - 8°C / 39°F - 46°F. Do not freeze! Bring to room temperature before use. Constant temperatures above 25°C / 77°F can reduce the shelf-life of the product.

12. Batch number and expiry date

The batch number should be specified to identify products in the case of enquiries. Products should no longer be used once the expiry date has elapsed.

13. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. For dental use only. Commercially available medical gloves do not provide protection against the sensitization effect of methacrylates. If the product comes into contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice.

14. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

15. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination.

16. Warranty

The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

17. Other notes for Europe

If the user and/or patient become aware of serious incidents connected with application the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides. The Summary of Safety and Clinical Performance of SAREMCO Products can be found in the European database on medical devices (EUDAMED-<https://ec.europa.eu/tools/eu-damed>).

18. Production / distribution

SAREMCO Dental AG
Gewerbestrasse 4
CH-9445 Rebstein / Switzerland
Tel: +41 (0) 71 775 80 90
Fax: +41 (0) 71 775 80 99
info@saremco.ch
www.saremco.ch

Edited 10-2022 | D600150

Class IIa medical devices
Class III medical devices (Canada)

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