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Instruction for use
cmf etch

- 1. Product description**
Dentine-preserving, buffered, phosphoric acid-containing etching gel for the total-etch technique.
- 2. Composition**
Water, phosphoric acid (35 % wt), phosphoric salt, gel former, colorant
- 3. Intended Use**
Sarecco dental bonding agents promote adhesion of direct or indirect prosthesis for reconstruction or correction of functionally compromised natural dentition (e.g., deficient teeth).

4. Indication
Acid etching of enamel and dentin surfaces prior to the application of primers, adhesives, or sealers, etc.

5. Contra-indication
Opened pulp, pulpitis, known allergy to methacrylates.

6. Patient target group
cmf etch is suitable for use in all patients without any age or gender restrictions.

7. User
cmf etch should only be used by a professionally trained dental practitioner.

8. Side effects
In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

9. Interactions
Underfilling materials based on glass ionomer and calcium hydroxide are etched by **cmf etch**.

10. Processing stages

10.1. Drying
Due to the following application of the filling material, a drainage of the working field is required. Avoid any contamination of the cavity with blood or saliva. A rubber dam is recommended.

10.2. Cavity preparation
Prepare the cavity as usual. Clean and dry the cavity.

10.3. Etching of tooth substance (total-etch technique)

10.3.1. Use a new syringe tip for each patient. Before using in the mouth, put a small amount of the gel on a cotton wool pad to familiarise yourself with the viscosity.

10.3.2. Apply **cmf etch** to enamel and dentine. Allow to have an effect on enamel and dentine for 30 seconds. The buffered etching gel reduces the risk of excessive demineralisation. Rinse and dry thoroughly. Avoid over-drying. Sufficiently etched enamel appears matt and white.

10.4. Adhesion
The adhesive has to be conditioned and applied in correspondence to the instruction for use of the product used. SAREMCO recommends **cmf adhesive system**, **els duobond**, **els unibond** and **apt unibond**.

10.5. Refilling the cmf etch syringe 2,5 ml
If necessary, **cmf etch** 2,5 ml can be refilled with **cmf etch** 25 ml. To do that, remove the caps from both syringes. Place the adapter supplied firmly on the 25 ml syringe and on the 2,5 ml syringe and screw tight. Once the syringes are firmly in place, slowly withdraw the plunger of the 2,5 ml syringe to collect the desired amount of material. Avoid unintentional leakage of the material. Carefully separate the syringes from the adapter. When not in use, always keep the syringes closed with the caps.

11. Storage
cmf etch has been developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Bring to room temperature before use.

12. Batch number and expiry date
The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

13. Precautionary measures
Close containers after each use with the right lid. **cmf etch** is corrosive! Avoid contact with skin, oral mucosa, eyes, clothes, metal etc. Keep out of reach of children. For dental use only. The principles of good personal and safe working practice apply.

14. Emergency measures
In case of direct contact with skin, oral mucosa or eyes, rinse immediately and thoroughly with running water. In case of contact with the eyes, consult an eye specialist.

15. Hygiene
Use application instruments for one patient

only. Dose products away from patients to avoid contamination.

16. Warranty
The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

17. Other notes for Europe
If the user and/or patient becomes aware of serious incidents in connection with the use of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides. The Summary of Safety and Clinical Performance of SAREMCO Products can be found in the European database on medical devices (EUDAMED-https://ec.europa.eu/tools/eudamed).

18. Production / distribution
SAREMCO Dental AG
Gewerbestrasse 4
CH-9445 Rebstein / Switzerland
Tel: +41 (0) 71 775 80 99
Fax: +41 (0) 71 775 80 99
info@sarecco.ch
www.sarecco.ch

Edited 10-2022 | D600165

Class IIa medical devices
Class III medical devices (Canada)

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2. Composition
Water, phosphoric acid (35 % wt), phosphoric salt, gel former, colorant

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