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2. Composition

inorganic fillers (barium glass and silica 76 % wt, 52 % by volume, particle size between 0.04 - 3.00 µm, median 0.7 µm), BisGMA, BisEMA, TEGDMA, catalysts, inhibitors, pigments

3. Intended Use

Saremo restoratives are intended for the direct reconstruction or correction of functionally compromised natural dentition (e.g., deficient teeth). Restorative material intended for the reconstruction or correction of natural dentition.

4. Indication

- Restoration of class I, II, III, IV and V cavities on anterior and posterior teeth
- Extra orally tempered inlays and onlays (direct and indirect), veneers
- Aesthetic corrections of interdental spaces, enamel hypoplasia, discolorations
- Blockings, fractures

5. Contra-indication

Opened pulp, pulpitis, known allergy to methacrylates.

6. Patient target group

mhc composite is suitable for use in all patients without any age or gender restrictions.

7. User

mhc composite should only be used by a professionally trained dental practitioner.

8. Side effects

In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

9. Interactions

Avoid underfilling materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

10. Processing stages

10.1. Tooth cleaning

Brush the tooth being treated and the teeth next to it with fluoride-free toothpaste. If necessary, clean interdental spaces using strips and dental floss.

10.2. Shade selection

SAREMCO shades relate to VITA standards.

10.3. Drying

Dry sufficiently, apply rubber dam.

10.4. Cavity preparation

Prepare the cavity as usual. Undercuts and beveled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity.

10.5. Underfilling

In case of near-to-pulp preparations, it is advisable to protect the pulp using an underfilling (e.g. calcium hydroxide). Do not use zinc oxide eugenol!

10.6. Adhesion

The adhesive has to be conditioned and applied in correspondence to instruction for use of the product. SAREMCO recommends *James-2, els unbond* and *cmf adhesive system*.

10.7. Restoration with mhc composite

After each usage, turn the syringe screw back by one turn and re-attach the screw cap. Gradually adapt *mhc composite* up to a maximum layer thickness of 2 mm. When using matrices, remove them and expose filling to light again from lingual and buccal.

Warning: Protect light-curing products from strong sources of light.

The specified exposure time refers to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm² and a wavelength of 400 - 500 nm. They apply to a maximal layer thickness of 2 mm. The necessary exposure time may vary depending on the light source and its instructions for use. In case of doubt, check the lamp's light output and the necessary exposure time before operation in vitro. During polymerization, an inhibition layer forms on the surface; this must not be touched or removed if other composite layers are going to be applied.

Exposure time for increments of 2 mm

light power	≥ 500 mW/cm ²	≥ 800 mW/cm ²
mhc opaque color	40 sec	30 sec
mhc all other colors	20 sec	20 sec

10.8. Finishing, polishing

Prepare the filling with 40 µm and 12 µm diamond burs. Polish to a high gloss using polishing brushes, polishing discs, strips or silicone polishers. The filling can be prepared and polished immediately after polymerization.

11. Storage

Protect light-curing products from strong sources of light and heat! *mhc composite* was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Bring to room temperature before use.

12. Batch number and expiry date

The batch number should be specified to identify products in the case of enquiries. Products should no longer be used once the expiry date has elapsed.

13. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. For dental use only. Commercially available medical gloves do not provide protection against the sensitization effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice. When polishing or removing composites, it is recommended to always use a water-cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

14. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

15. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination.

16. Warranty

The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

17. Other notes for Europe

If the user and/or patient becomes aware of serious incidents in connection with the use of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides. The Summary of Safety and Clinical Performance of SAREMCO Products can be found in the European database on medical devices (EUDAMED-<https://ec.europa.eu/tools/eudamed/>).

18. Production / distribution

SAREMCO Dental AG
Gewerbstrasse 4
CH-9445 Rebstein / Switzerland
Tel: +41 (0) 71 775 80 90
Fax: +41 (0) 71 775 80 99
info@saremco.ch
www.saremco.ch

Edited 11-2022 | D600172

Class IIa medical devices
Class III medical devices (Canada)

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