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saremco print - DENTURETEC

1. Product description

DENTURETEC is a light-curing, flowable polymer based with methacrylic acid ester for production of 3D-printed prosthesis bases.

2. Composition

BiEMA, Urethanmethacrylate, TEGDMA, pyrogenic silica, catalysts, inhibitors, pigments.

3. Intended use

saremco print products provide light-curing 3D printable resin-based materials for the correction or reconstruction of functionally compromised natural dentition (e.g., missing teeth or deficient teeth) by manufacturing of customized 3D-printed dental prostheses.

4. Indication

Production of prosthesis bases with the aid of the 3D-printing machine «ASIGA MAX UV».

5. Contra-indication

Do not use the product in case of a known allergy to one or more ingredients.

In case of doubt, clarify and exclude a possible allergy with the help of a specific allergy test before using **DENTURETEC**.

DENTURETEC must not be used for any other purposes than prosthesis bases. Any deviation from this instruction for use may have negative effects on the chemical and physical quality of the prosthesis bases produced from **DENTURETEC**.

6. Patient target group

DENTURETEC is suitable for use in all patients without any age or gender restrictions.

7. User

DENTURETEC should only be used by a professionally trained dental practitioner.

8. Side effects

In individual cases, contact allergies to products with similar composition have been described.

9. Interactions

None known.

10. Processing stages

10.1. Generating printing file
Generate the printing file of the desired prosthesis base by using appropriate software (Composer) and deliver it suitably to the printer. Please observe the corresponding instruction for use of software and printer. It is recommended to set up the prosthesis base vertically. Make sure that enough supports are generated

10.2. Printing

Work as clean as possible, as dirty reservoirs or machines can cause deformation/discolouration and therefore failure of the printed objects. Briefly shake the liquid material and pour it into the reservoir of the 3D-printing machine. The corresponding printing parameters are stored in the data base on the ASIGA website and must be loaded before printing. Start the printing process by following the instruction for use of the printer. A working temperature of 35°C / 95°F is required. It is always recommended to wait until the printer has reached this temperature.

10.3. Cleaning

After the printing process is completed, remove the building platform from the machine. During removing the prosthesis base and the following cleaning steps, wearing gloves (nitrile gloves) and protective goggles are advised.

Place the platform on a piece of paper or cloth with the built jobs facing upwards. Remove the printed jobs from the platform by using a suitable instrument (putty knife). To remove excess material, clean the printed job with an alcohol-soaked (96%) cloth and possibly a brush soaked in an alcohol solution until all resin remains are completely removed. Then dry the prosthesis base thoroughly with an air syringe.

Warning: Protect light-curing products from strong light sources!

10.4. Finishing the prosthesis base

To achieve the desired material properties and biocompatibility, post-curing of the completely dried and cleaned printed objects is necessary. For final polymerisation place the printed jobs in a UV-light box.

Note: time of curing depends greatly on type of lamps / lightbox used. The final properties and the final colour depend on the post-curing process. Post-curing is a UV-light treatment to ensure that **saremco print** materials obtained full polymer conversion, the residual monomer is reduced to a minimum and the highest mechanical properties are achieved. This procedure is a necessary step to attain a biocompatible end-product.

It is suggested to use the polymerisation unit "Signum HiLite Power" from Heraeus Kulzer (2 x 180s) or the UV-Flash device "Otoflash G171" from NK-Optik (4000 flashes), the latter device being preferred due to curing under nitrogen. It is recommended to let the material cool down a little between the curing cycles. In general, all lightboxes for light-curing veneering materials can be used that cover a wavelength range of 320 - 500 nm.

Lightboxes with integrated flash light allow shorter exposure time compared to conventional lamps. Always follow the respective instruction of use of the polymerisation unit. Blast the surface of the prosthesis base with blast polishing material carefully (e.g. Perlablast micro BEGO). Afterwards remove the support structures by using a cut-off-wheel or a cutter. Wearing nitrile

gloves, safety glasses and dust mask is advised during this finishing process.

10.5. Inserting of artificial teeth in the printed, prefabricated prosthesis base

Roughen the base surface of the printed artificial teeth for example by sandblasting (Al2O3, 110 µm), apply a primer and a fixing material, insert in the prosthesis according to the natural shape and polymerise. Alternatively, **DENTURETEC** can also be used as fixing material. For this purpose, put a small amount of material with a brush on the roughened teeth-surface of the artificial tooth, put it into the prosthesis, eliminate any excess material and light cure it from all sides for at least 20 seconds. The polymerisation lamp should work in the range of 320 to 500 nm and have at least a light output of 600 mW/cm².

11. Storage

Protect this product from strong light and heat sources! The recommended storage temperature is between 4°C and 28°C / 39°F and 82°F. Close the package after each use.

12. Batch number and expiry date

The batch number is used to identify the product in case of queries. Do not use this product after the expiration date.

13. Precautionary measures

For dental use only. Keep out of reach of children. The use of nitrile gloves while working with **DENTURETEC** is recommended until post-curing. Commercially available medical gloves do not provide effective protection against the sensitising effect of methacrylates. If the product comes into contact with the glove, remove the glove and dispose it of, wash your hands immediately with soap and water and put on a new glove. In case of an allergic reaction, consult a doctor. When polishing or removing composites, it is recommended to always use a water cooling system and a good extraction system, to ventilate the dental laboratory frequently and to wear masks with high particle filtration efficiency for small particle sizes.

14. Emergency measures

In case of direct contact of the uncured material with the oral mucosa, rinse with water. In case of contact with eyes, rinse thoroughly with water and consult an eye specialist.

15. Hygiene

Prosthesis bases made of **DENTURETEC** should not be cleaned with chemical products. Cleaning with water is sufficient. The finished restorations can - if necessary - be disinfected with an ethanolic solution.

16. Warranty

Our liability is restricted to the quality of our products. In case of a product being of defective quality only its value is replaced. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental.

The Summary of Safety and Clinical Performance of SAREMCO Products can be found in the European database on medical devices (EUDAMED-https://ec.europa.eu/tools/eudamed).

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