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1. Produktbeschreibung

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Instruction for use peek primer

1. Product description

peek primer is a light-curing, low-viscosity bonding promoter for artificial teeth and veneers based on polymethyl methacrylate (PMMA), high-performance polymers such as polyaryletherketone (PAEK) or polyethyletherketone (PEEK) and methacrylate-based composites. **peek primer** ensures a firm bond when bonding with (meth)acrylate-based materials and is designed for use in the dental laboratory for indirectly manufactured dentures.

2. Composition

Triallyl cyanurate, Bis-GMA, pentaerythritol triacrylate, initiators, catalysts, inhibitors

3. Intended Use

Saremco cements and auxiliary bonding promoters are intended for the indirect reconstruction or correction of functionally compromised natural dentition (e.g., deficient teeth).

4. Indication

- Priming of artificial teeth and veneers made of PMMA or composite base
- Priming of frameworks, prosthesis bases as well as crowns and bridge materials based on PMMA, composite base and thermoplastic high-performance polymers such as PAEK, PEEK and similar.
- Repair of crowns and bridges made of plastic respective composite.

5. Contra-indication

Known allergy to methacrylates. **peek primer** is not suitable for conditioning metal or ceramic or rather zirconium dioxide frameworks.

6. Patient target group

peek primer is suitable for use in all patients without any age or gender restrictions.

7. User

peek primer should only be used by a professionally trained dental practitioner.

8. Side effects

In individual cases, contact allergies have been described to products of a similar composition.

9. Interactions

Avoid contact with materials which may hinder polymerization owing to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

10. Processing stages

10.1. Conditioning

Blast the materials to be processed with 110 µm aluminum oxide at 2 - 3 bar (roughening, increasing the surface) and clean with oil-free compressed air. Then use a brush to apply the **peek primer** to the surface to be conditioned and polymerize in the light-curing device according to point 10.2. As an alternative to blasting, the surfaces can also be roughened using abrasives. Do not steam jet the object. Remove dirt with alcohol and a clean brush.

Note: Apply **peek primer** thinly, apply only once and polymerize immediately.

10.2. Light polymerization devices and polymerization times

Different light sources are used in light polymerization devices (laboratory & practice). The necessary exposure time can therefore vary depending on the light source and its instructions for use. All light sources with a wavelength range of 430 to 470 nm are suitable. In case of doubt, check the light output of the lamp and the necessary exposure time before treatment in vitro.

Halogen-lamps with 700 mW/cm ²	40 sec
LED lamps with 1100 mW/cm ²	20 sec
HiLight power	90 sec
Otoflash G171 from NK-Optik	2'000 flashes

Note: Complete polymerization is guaranteed when the surface is completely dry after curing.

11. Storage

Protect light-curing products from strong sources of light and heat! **peek primer** was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Constant temperatures above 28°C / 82°F can reduce the shelf life of the product.

12. Batch number and expiry date

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

13. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. For dental use only. Commercially available medical gloves do not provide protection against the sensitization effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice.

Hazard warnings: May cause an allergic skin reaction. Causes serious eye damage. Harmful to aquatic organisms, with long-term effect.

14. Emergency measures

In case of contact with the eyes, rinse cautiously with water for a few minutes. Remove any contact lenses, if possible. Continue rinsing. Consult an eye specialist.

15. Hygiene

Use application instruments for one patient only. Dose products away from patients to avoid contamination. The usual hygiene rules

must be observed. Dispose of the content/ container in accordance with local, regional, national, and international regulations.

16. Warranty

The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

17. Other notes for Europe

If the user and/or patient becomes aware of serious incidents in connection with the use of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides.

The Summary of Safety and Clinical Performance of SAREMCO Products can be found in the European database on medical devices (EUDAMED-<https://ec.europa.eu/tools/eudamed/>).

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