

# CONTEC



Free of  
BisGMA\*  
TEGDMA\*  
HEMA\*!



The secure connection.

# The light-curing bracket adhesive.

For the adhesion of metal, ceramic and plastic brackets and buccal tubes.

CONTEC Ic offers all the advantages of a light-curing adhesive. This gives you more time to place the brackets before the adhesive cures.

The high viscosity ensures that the adhesive stays on the bracket and does not drift. At the same time, the bond strength of CONTEC Ic is excellent. It is therefore the solution for all bonding tasks.

The ergonomically-formed syringe makes for easy handling and good control over the amount of adhesive used.

The special advantage of CONTEC Ic is that it is free of BisGMA\*, TEGDMA and HEMA. These materials have been proven to have allergenic potential.



#### CONTEC Ic set (REF 163-110-00) consisting of:

- 3 x 3 g CONTEC Ic adhesive
- 1 x 5 ml CONTEC Ic/lcr primer
- 2 x 2.5 ml CONTEC etching gel
- 6 x CONTEC disposable tubes for etching gel
- 50 x CONTEC disposable brushes
- 1 x CONTEC combination instrument
- 1 x CONTEC mixing pad

For approx. 800 brackets (40 treatment cases).



\* When used with metal and ceramic brackets, CONTEC Ic is free of BisGMA (bisphenol A glycerolate dimethacrylate), TEGDMA (triethylene glycol dimethacrylate) and HEMA (hydroxyethyl methacrylate). CONTEC Ic, used with plastic brackets, is free of TEGDMA and HEMA.

- Excellent bond strength
- Easy-to-dose adhesive paste with very good stability
- Cures under all bracket bases (an additional primer is needed for plastic brackets)
- Free of BisGMA\*, TEGDMA\* and HEMA\*



Single pack	Contents	REF
CONTEC Ic adhesive	2 x 3 g syringes	163-111-00
CONTEC Ic/Icr primer	1 x 5 ml bottle	163-113-10
CONTEC plastic primer	1 x 5 ml bottle	163-100-10
CONTEC etching gel	2 x 2.5 ml syringes	163-115-00
CONTEC disposable tubes for etching gel	10 pieces	163-119-00
CONTEC disposable brushes	50 pieces	163-117-00
CONTEC combination instruments	5 pieces	163-118-00
CONTEC mixing pads	2 pieces	163-116-00

Developed in cooperation with:

**SAREMCO**  
DENTAL

# The light-curing retainer adhesive.

For the adhesion of fixed retainers.



CONTEC lcr has been developed to fulfill the special requirements for bonding fixed orthodontic retainers. With CONTEC lcr, retainers can be secured safely and with the maximum patient comfort possible. Yet the adhesive is durable and lasting and has excellent bonding characteristics.

The slightly bent tubes and the delicate syringes make handling in the lingual area of the anterior teeth easy.

CONTEC lcr is free of TEGDMA\* and HEMA\*. CONTEC lcr is therefore the first choice for patients who have been diagnosed with an allergy against one of these substances.

Single pack	Contents	REF
CONTEC lcr adhesive	1 x 3,4 g syringe	163-131-00
CONTEC lcr disposable tubes	10 pieces	163-120-00
CONTEC lcr/lcr primer	1 x 5 ml bottle	163-113-10
CONTEC etching gel	2 x 2.5 ml syringes	163-115-00
CONTEC disposable tubes for etching gel	10 pieces	163-119-00
CONTEC disposable brushes	50 pieces	163-117-00
CONTEC combination instruments	5 pieces	163-118-00
CONTEC mixing pads	2 pieces	163-116-00

Developed in cooperation with:

**SAREMCO**  
DENTAL

\* CONTEC lcr is free of TEGDMA (triethylene glycol dimethacrylate) and HEMA (hydroxyethyl methacrylate).

- Pasty adhesive, high filler content
- Exceptionally durable and lasting
- Cures to form an ultra-smooth surface
- Free of TEGDMA\* and HEMA\*

Free of  
TEGDMA\*  
HEMA\*!




**CONTEC Icr set (REF 163-130-00) consisting of:**

- 2 x 3.4 g CONTEC Icr adhesive
- 6 x CONTEC Icr disposable tubes
- 1 x 5 ml CONTEC Icr/Icr primer
- 2 x 2.5 ml CONTEC etching gel
- 6 x CONTEC disposable tubes for etching gel
- 50 x CONTEC disposable brushes
- 1 x CONTEC combination instrument
- 1 x CONTEC mixing pad

For approx. 80 retainers, twisted or approx. 250 retainers 3-3.

# The ready-to-use, chemically-curing adhesive.

For adhering metal, ceramic and plastic brackets and buccal tubes.



ConTec Go! is a chemically curing, ready-to-use bracket adhesive consisting of an adhesive paste and an activator. The curing process starts by means of pressure. This allows the brackets to be positioned on the tooth before the adhesive cures.

There is sufficient strength in the bond after only 5 minutes for the archwire to be ligated.

ConTec Go! has a high viscosity and a consistency that guarantees optimum stability on the bracket. ConTec Go! also maintains its color throughout the period of treatment.

Single pack	Contents	REF
ConTec Go! adhesive	2 x 2.5 g syringes	163-321-00
ConTec Go! activator	1 x 10 ml bottle	163-314-00
ConTec Go! etching gel	2 x 2.5 g syringes	163-312-00
ConTec Go! disposable tubes for etching gel	10 pieces	163-322-00
ConTec Go! disposable brushes	50 pieces	163-317-00
ConTec Go! brush holders	5 pieces	163-319-00
ConTec Go! spatulas	50 pieces	163-318-00
ConTec Go! mixing pads	2 pieces	163-316-00

- Viscous adhesive for markedly higher stability
- Curing time reduced to 5 minutes
- Ergonomically shaped syringes for comfortable handling and safe dosing

**ConTec Go! Set (REF 163-307-00) consisting of:**

- 3 x 2.5 g ConTec Go! adhesive
- 1 x 10 ml ConTec Go! activator
- 2 x 2.5 g ConTec Go! etching gel
- 6 x ConTec Go! disposable tubes for etching gel
- 50 x ConTec Go! disposable brushes
- 1 x ConTec Go! brush holder
- 1 x ConTec Go! spatula
- 1 x ConTec Go! mixing pad

For approx. 600 brackets (30 treatment cases).



# Dentaurum

Germany | Benelux | España | France | Italia | Switzerland | Australia | Canada | USA  
and in more than 130 countries worldwide.



DENTAURUM  
QUALITY  
WORLDWIDE  
UNIQUE

	<a href="http://www.dentaurum.com">www.dentaurum.com</a>		

Date of information: 10/19  
Subject to modifications

