

THE PERMANENT CROWN FROM THE 3D PRINTER

saremco print - CROWNTEC

THE „GAME CHANGER“

Teeth and denture base
customized from the
3D-printer with
CROWNTEC and
DENTURETEC*

CROWNTEC is a qualitatively outstanding composite resin for the manufacturing of highly biocompatible permanent restorations.

- for permanent crowns, inlays, onlays, veneers and artificial teeth
- highly esthetic one-phase material with perfectly matched opacity
- outstanding biocompatibility: no elutable substances in an aqueous environment
- ready to use without requiring a lot of shaking
- refining possible with saremco els paintart



*see separate product information

More information: www.saremco.ch

Technical Data

Consistency	middle viscous, homogeneous
Color	according to VITA® standard
Depth of cure (DIN EN ISO 4049)	≥ 1,5 mm
Flexural strength (DIN EN ISO 4049)	> 130 MPa
E-Modulus (DIN EN ISO 4049)	> 4000 MPa
Barcol hardness	> 40
Shelf life	3 years
Medical device	class IIa

NEW PRINT- FILE
FOR ASIGA!
NOTE THE STICKER

Order information (500g bottle)

REF	Product Description	Availability
8063	CROWNTEC, A1	on stock
8052	CROWNTEC, A2	on stock
8051	CROWNTEC, A3	on stock
8055	CROWNTEC, A3.5	on stock
8057	CROWNTEC, A4	on demand
8065	CROWNTEC, B1 New!	on stock
8056	CROWNTEC, B2	on demand
8053	CROWNTEC, B3	on demand
8058	CROWNTEC, C3	on demand
8059	CROWNTEC, C4	on demand
8054	CROWNTEC, D3	on demand
8066	CROWNTEC, SW (bleach-color) New!	on stock

