

# saremco print | Compatibility Overview

| 3D PRINT MATERIAL   | 3D PRINTER  | CLEANING  | POST-TREATMENT   |
|---|---|---|--|
| <p><b>DENTURETEC</b><br/>for all types of denture bases</p> | <p><b>ASIGA</b><br/>all 385nm printers</p> <p>More devices follow</p> | <ul style="list-style-type: none"> <li>• Cleaning by hand (brushes and cloths) with minimal use of IPA (96%)</li> <li>• Air dry the surface and inside of the print object</li> </ul> | <ul style="list-style-type: none"> <li>• <b>OtoFlash</b> (NK-Optik)<br/>2 x 2000 Flashes</li> <li>• <b>HiLite Power</b> (Kulzer)<br/>2 x 180 s</li> <li>• More devices follow</li> </ul> |

Please note: This compatibility overview does not replace the instructions for use. Please read the instructions for use carefully.

Edited: 04/2022 | D600245 | Class IIa Medical Device

**More information at [saremco.ch/en/downloads/instruction-for-use/](https://saremco.ch/en/downloads/instruction-for-use/)**